

Mandatory Internship Report



Hochschule Rhein-Waal

and



Siemens Healthineers Innovation Think Tank Mechatronic Products

Faculty of Technology and Bionics

Mechatronics Systems Engineering

16th August ,2024

Syed Mustafa Mohiuddin (26293)

Supervised by:

Dr. Kerim Kara, Siemens Healthineers, ITT Technology Implementation Lead

Dr.-Ing. Stéphane Danjou, Hochschule Rhein-Waal

Table of Contents

Confidentiality Clause.....	4
Declaration of Authenticity.....	5
Internship Report Approval Page.....	6
1 Foreword and Acknowledgements.....	7
2 Company	8
2.1 Profile	8
2.2 Innovation Think Tank at Siemens Healthineers.....	10
2.3 Products and Market	11
3 Internship Description	12
3.1 Activities	12
3.2 Working Conditions and Functions	12
4 PROJECTS.....	13
4.1 Active Collision Avoidance System for Nexaris Angio-CT Sliding System	13
4.1.1 Project Description.....	13
4.1.2 Task Description	14
4.1.3 Tools / Methods Used	14
4.1.4 Project Goals	15
4.1.5 Project Approach.....	16
4.1.6 Testing stage:	20
4.1.7 Conclusion.....	23
4.2 Collimator Rotation Angulation Detection	24
4.2.1 Project Description.....	24
4.2.2 Task Description	25
4.2.3 Tools / Methods Used	28
4.2.4 Project Goals	29
4.2.5 Project Approach.....	30
4.2.6 Testing Stage.....	31

4.2.6 Conclusion.....	31
4.2.7 Key Achievements include:	32
5 Professional Insights	33
5.1 Technical and Practical Knowledge.....	33
5.2 Non-technical outcomes.....	33
5.3 Conclusion.....	34
6 Certificates at ITT MP.....	35
7 Internship Certificate.....	36
8 Literature/ References	37

Confidentiality Clause

This internship report contains confidential data from Siemens Healthineers.

This work may only be made available to the first and second examiner and authorized members of the board of examiners at Hochschule Rhein-Waal. Any publications and duplication are prohibited.

An inspection of this work by third parties requires the expressed permissions of the authors and Siemens Healthineers.

Declaration of Authenticity

Declaration

I, Syed Mustafa Mohiuddin, declare that the research work presented here is from the best of my knowledge and belief, original and the result of my own investigations. The cooperation I got for this research work is clearly acknowledged. To the best of my knowledge, it does not contain any materials those are written by others or published already except mentioned with due references in the text as well as with the quotation marks. This work has not been published, submitted, either in part or whole intended for reward, degree at this or any other University.

Date, Location

Signature

Internship Report Approval Page

I, the undersigned, confirm that I have reviewed and evaluated the internship report submitted by Syed Mustafa Mohiuddin for the completion of Internship at Siemens Healthcare GmbH. After a thorough examination of the report's content and quality, I am satisfied with the work presented. I acknowledge that the report accurately represents the work completed during the internship and fulfills the academic or professional requirements set forth for this internship program.

Supervisor's Signature _____

1 Foreword and Acknowledgements

This internship report was created describing the tasks performed during my internship at Innovation Think Tank Mechatronic Products (ITT MP) located at Siemens Healthineers (SHS), Kemnath. The internship lasted from April 2024 to August 2024 as a part of the mandatory internship requirement by Hochschule Rhein-Waal.

I would like to express my sincere gratitude to Siemens Company for granting me the opportunity to work as an intern. The knowledge and skills I gained during my time there, including working with new technologies, creating physical prototypes, and managing projects, have been incredibly valuable. I want to extend a special thanks to Prof. Haider and Dr. Kerim Kara for giving me the chance to be part of their diverse team. I also want to acknowledge the constant support and assistance from the coordinators - Mr. Konstantinos Gkekas and Mr. Bernd Niklas Axer. Their willingness to share their knowledge and insights played a significant role in my professional growth. This internship at Siemens has given me real-world exposure to projects and has contributed greatly to my development as a professional. Moreover, I would also like to thank Dr.-Ing. Stéphane Danjou from Hochschule Rhein-Waal for being my supervisor and for providing me guidance and support at times of need, in addition I would like to thank my colleagues at ITT MP for such a warm welcome into the team and assisting me in projects.

2 Company

2.1 Profile

Siemens Healthineers, a pioneering global medical technology firm headquartered in Erlangen, Germany, is recognized globally for its transformative contributions to healthcare. With over 52,000 employees worldwide, the company is an industry leader, generating 14.5 billion Euros in revenue and retaining a profit of 2.5 billion Euros in 2019. As a company with more than 120 years of experience and 18,500 patents globally, SHS has nearly 66,000 employees in more than 70 countries that are shown in Fig. 1.1 to innovate and shape the future of healthcare [1]. Its strategic accomplishments underscore innovation and operational excellence. Siemens Healthineer's advanced medical imaging and diagnostic solutions make it an influential player in global healthcare. The company's diverse facilities are strategically distributed across continents, symbolizing its commitment to serving diverse markets. Through its unwavering dedication to innovation, Siemens Healthineers continues to redefine medical technology's impact on patient care and diagnostics.




Figure 1.1: Geographical Distribution of Siemens Healthineers Facilities [1]

The plant facilities located in Kemnath which has now become the world's largest location for mechatronics in the field of imaging medical technology has dated back to early 1962. The site has a complete process and competence chain, beginning in innovation, development, production, assembly, and system testing to commissioning in the hospital. With the guiding principle of the Kemnath location, "Emotion in motion", the site was awarded with "Factory of the Year" in 2012 for "Outstanding Change Management" for successfully restructuring of technology production with close involvement of the employees.



Figure 1.2: Siemens Healthineers Location in Kemnath.

2.2 Innovation Think Tank at Siemens Healthineers

Established in 2005, Innovation Think Tank (ITT) is part of the Chief Technology Office of SHS. ITT is an interdisciplinary and self-sustaining infrastructure within the organization that focuses on innovation by implementing the ITT Framework and is accessible for the company founded and lead by Professor Sultan Haider. As a global infrastructure, multiple innovation labs and Innovation Think Tank programs can be found across the globe addressing projects from within SHS and host institutions. ITT MP is one of the innovation labs located in the SHS facility in Kemnath. Comprised of a diverse and transdisciplinary team of designers and engineers from different fields, ITT MP has continuously contributed to defining, shaping, and the positing of the Product Management (PRM), Research and Development (R&D) and Supply Chain Management (SCM) through over 960 project modules with Siemens Healthineers Mechatronic Products (SHS MP) since 2010.

ITT aims to drive healthcare innovation to improve patient outcomes and experiences. By focusing on patient-centric solutions, ITT advances technologies that help healthcare providers deliver superior care. Their services include technology analysis, competitor insights, market trends, and solution development. ITT's skilled team uses data-driven methods to create transformative solutions in medical imaging, diagnostics, healthcare IT, and other areas.




Figure 1.3: Organizational structure of Innovation Think Tank in Siemens Healthineers [2]

2.3 Products and Market

Present in more than 70 countries around the globe, SHS portfolio is at the center of treatment pathways and clinical decision-making processes. The modalities provided are within the scope of

- **Imaging:** As a market leader in diagnostic imaging, SHS provides systems used for angiography, computed tomography (CT), fluoroscopy, magnetic resonance imaging (MRI), molecular imaging, x-ray, ultrasound and imaging software.
- **Laboratory and point of care Diagnosis:** SHS offers services that support the needs of laboratories and healthcare facilities in laboratory diagnostics in conjunction to automation, informatics, and services. The spectrum includes immunoassay, chemistry, hematology, molecular, and urinalysis testing.
- **Advanced Therapies:** Supporting advanced therapies, SHS provides high-end angiography systems, mobile C-arms and hybrid operating rooms to empower innovative therapy concepts and minimally invasive procedures.
- **Services:** In ensuring well performing systems, trained staff and optimized processes, SHS is a trusted partner to provide the customer's confidence in delivering the right diagnosis to patients.

3 Internship Description

3.1 Activities

During my internship at ITT MP, I was tasked with developing proof-of-concept solutions for healthcare challenges using ITT methodologies. My role involved contributing to multiple projects, including collision avoidance systems and angular detection for x ray-collimator. The internship provided a dynamic work environment, alternating between office-based tasks and remote work, facilitating a balance between collaboration and independent research.

3.2 Working Conditions and Functions

- Location: Siemens Healthineers, Kemnath, Germany
- Duration: 01. April.2024 – 16. August.2024 (20 Weeks)
- Working hours: 8:00 AM to 4:00 PM

The initial 2 weeks were dedicated to onboarding tasks and watching informative Learn4u videos. These videos offered detailed overview of medical procedures and Siemens-related offerings, providing a very productive understanding of the field. The work atmosphere was conducive to productivity and balanced work-life integration. Although some days demanded extra hours, they contributed to my overall enriching experience. Adjusting to this new work scenario was seamless, largely thanks to the supportive and welcoming team. Overcoming language and cultural differences was easier than anticipated, as my supervisors and coordinators readily provided guidance.

My internship encompassed a diverse range of impactful tasks, contributing significantly to my professional growth. Notably, I developed essential teamwork and communication skills. Working alongside colleagues from diverse backgrounds enhanced my collaboration abilities and effective communication across cultures.

This internship broadened my horizons and gave me practical skills for a global work environment. Following section goes into the projects.

4 Projects

4.1 Active Collision Avoidance System for Nexaris Angio-CT Sliding System

4.1.1 Project Description

The Active Collision Avoidance System (ACAS) for sliding CT gantries represents a significant advancement in medical imaging safety and efficiency. Sliding CT gantries, essential for modern diagnostic and interventional procedures such as CT scans and angiography, require seamless mobility. However, their dynamic movement introduces safety risks including potential collisions with equipment, healthcare personnel, or patients. These risks can lead to equipment damage, increased maintenance costs, workflow disruptions, and safety concerns.



Figure 1.4: Nexaris Angio-CT Sliding System [3]

The ACAS project was initiated to address these challenges by developing an advanced sensor-based system capable of real-time obstacle detection and response. The goal was to enhance safety, improve operational efficiency, and reduce maintenance costs while providing a user-friendly solution that integrates seamlessly with existing CT gantry systems.

To overcome these limitations, the project aimed to implement an advanced sensor-based collision avoidance system capable of real-time obstacle detection and response. By integrating cutting-edge sensors and intelligent software, the project sought to enhance the safety, efficiency, and reliability of sliding CT gantry operations, providing a more secure and streamlined experience for medical professionals and patients alike.

4.1.2 Task Description

The project involved several key tasks. The project was already in development and a new collision avoidance feature as to be added. Initially analyze the medical imaging workflow environment to identify potential collision scenarios and safety risks, considering the movement of equipment and personnel. Next, develop detailed use cases to guide the design of the collision avoidance system, covering various room layouts, obstacles, and operator interactions. Select and configure a suitable sensor, for real-time obstacle detection, and integrate it with the sliding CT gantry. Design a 3D model of the sensor mounting bracket using Siemens NX CAD software and 3D printing the designed part, including designing an electrical schematic and soldering for connecting the sensor to the indicator. Integrate the system, conduct extensive testing to ensure accurate obstacle detection, and minimize false alarms. Finally, deploy the system in real-world medical imaging settings, gather operator feedback, and refine the system to enhance its performance.

4.1.3 Tools / Methods Used

1. Sensors:
 - Nano Scan 3 Pro from SICK: 2D Lidar Sensor chosen for its accuracy and reliability in real-time obstacle detection for our developed use cases.
2. Configuration Software:
 - Safety Designer Software: Used for configuring the detection field and the electrical schematic and ensuring its effective integration with the CT gantry.

- Siemens NX: Utilized for designing the 3D model of the sensor mounting bracket.
3. Tools:
 - Soldering wires
 - Ultimaker Cura for 3D Printing
 4. Testing Methods:
 - Simulations: Tested the system under various scenarios to assess performance and reliability.
 - Field Testing: Deployed the system in real-world environments to validate its functionality and gather operator feedback.

4.1.4 Project Goals

The ACAS project was designed with several key goals in mind, all focused on improving safety, efficiency, and reliability in medical imaging environments. The primary goal was to enhance safety by reducing the risk of collisions between the sliding CT gantry and the possible obstacles in the Operations Room, thereby ensuring a safe environment for patients and staff. This required the implementation of a system capable of detecting and responding to obstacles in real-time, effectively minimizing accidents and injuries.

Improving efficiency was another critical goal, aimed at minimizing workflow disruptions caused by false alarms or abrupt gantry stoppages. To achieve this, the system was optimized to accurately distinguish between true obstacles and benign objects, allowing for smooth and uninterrupted operation.

Reducing maintenance costs was also a key objective. By implementing collision avoidance measures that prevent mechanical stress and wear, the project sought to protect equipment from damage and extend its operational lifespan, ultimately reducing repair needs.

4.1.5 Project Approach

1. Initial Analysis and Design:

- Conduct a thorough workflow analysis to identify potential collision scenarios and define system requirements.
- Develop use cases to guide system design and functionality.

2. Sensor Selection and Configuration:

- Researched and evaluated sensor technologies and selected the Nano Scan 3 Pro.
- Configure the sensor using Safety Designer software and integrate it with the CT gantry.

3. Design and Implementation Stack light of Mounting Bracket:

- Design and print a 3D model for the sensor mounting bracket using Siemens NX and Ultimaker Cura for 3d printing , ensuring flexibility and precision in sensor positioning.

4. Integration and Testing:

- Soldering the sensor system with an indication stack light and performing simulations to test performance.
- Conduct field testing in medical imaging environments to validate effectiveness and gather operator feedback.

5. Refinement and Finalisation:

- Refine the system based on field test results and feedback.
- Finalize the design and prepare for deployment.




Figure 1.5: Sensor Selection Grid.

Possible Scenarios						
Detection Field		Ceiling mounted equipment : - Ceiling monitors - Ceiling OT lights - Radiation shield	Floor mounted equipment : - Roller container - IV stand - Anesthesia Equipment - Radiation shield	Moving personnel: - Gantry operator - Patient - Medical staff	Lying / Fallen objects : - Gloves - Dust Bins - Tissues - Wires - Components in OR(sharp instruments,head rest)	Reverse Movement of going back to rotation area
	3D Lidar #1 (with Brush)	✓	✓	✓	✓	✗
	3D Lidar #2 (with Brush)	✓	✓	✗	✓	✓
	3D Lidar #3 (with Brush)	✓	✓	✓	✓	✗
	2D Lidar (with Brush)	✓	✓	✓	✓	✗

Figure 1.6: Sensor Selection Research.

The sensor I selected was nano Scan 3 pro from SICK. Once the sensor was selected, it was configured using Safety Designer software and integrated seamlessly with the sliding CT gantry. This phase required meticulous calibration to ensure optimal sensor coverage and accuracy. Integration with existing gantry control systems was also performed to allow the sensors to communicate effectively with the gantry’s operational software.




Figure 1.7: nano Scan 3 pro [4]

Extensive sensor testing followed, verifying the sensor's performance in accurately and reliably detecting obstacles. The tests included simulations of various scenarios outlined in the use cases, ensuring that the system could handle different obstacles and conditions without false alarms or missed detections .



Figure 1.8: Safety Designer Software logo.[5]




Figure 1.9: Detection field configuration on the gantry.




Figure 2.0: Operation Room layout with ACAS.

4.1.6 Testing stage:




Figure 2.1: nanoScan3 configuration with stack light



Figure 2.2: Testing Setup with an obstacle and indication stack lights.

A crucial part of this integration involved the design and development of a 3D model for a mounting bracket for the sensors. Using Siemens NX, a sophisticated CAD software, then print a 3D model that allows the bracket to support three degrees of freedom (3 DOF). This design enables precise alignment and positioning of the sensors, ensuring they can capture the necessary data from various angles. The flexibility of the 3D model facilitates easier installation and adjustments in the field, accommodating different gantry configurations and enhancing the system's adaptability.



Figure 2.3: Siemens NX Logo.





Figure 2.4: Mounting Bracket CAD model.





Figure 2.5: Second Mounting Bracket CAD model.

The final phase of the project involved testing the system in real-world environments, deploying the system in actual medical imaging settings to validate its performance under real-world conditions. This testing phase provided valuable insights into the system's effectiveness and usability, allowing for further refinements and adjustments based on feedback from operators and stakeholders.

4.1.7 Conclusion

The Active Collision Avoidance System for Nexaris Angio- CT Sliding project successfully addressed critical safety and efficiency concerns in medical imaging environments. The project achieved its goals by enhancing safety, reducing maintenance costs, and improving workflow efficiency through advanced sensor technology and intelligent system design. The structured approach, including thorough analysis, precise sensor selection, and extensive testing, ensured that the ACAS met the demands of modern medical facilities. The integration of the Nano Scan 3 Pro sensor and the design of a customizable mounting bracket demonstrated the potential of advanced technologies in transforming healthcare environments.

Overall, the ACAS project exemplifies the effectiveness of combining cutting-edge sensor solutions with thoughtful system design to enhance safety and operational efficiency in critical medical imaging applications.

4.2 Collimator Rotation Angulation Detection

4.2.1 Project Description

The primary objective of the project was to develop a system for detecting and visualizing the rotation angulation of a x-ray collimator. Collimators are essential in various fields, particularly in medical imaging and radiation therapy, where precise alignment is crucial for accurate results. For this project, the Bosch Cross Domain Development Kit 110 (XDK) sensor kit was employed to monitor the orientation of the collimator. The goal was to create a reliable system that not only measures the rotation but also integrates multiple sensors to enhance accuracy and visualization of the data.




Figure 2.9: x-ray collimator movement. [6]



Figure 3.0: x-ray collimator setup. [6]

4.2.2 Task Description

The project was divided into several key phases:

1. Understanding the Infrastructure and Sensor Familiarization:

- Understand the concept of why we need and angulation detection for x-ray collimator and develop its use cases.
- Acquired detailed knowledge of the Bosch XDK110 sensor kit, including its technical specifications and operational principles.
- Set up the sensor and familiarized myself with its software environment and communication protocols.




Figure 3.1: XDK 110 capabilities. [7]

2. Code Development for Orientation Visualization:

- Developed code to process sensor data and visualize the orientation of the collimator using a combination of Gyroscope, Accelerometer and Magnetometer. This involved interpreting raw data from the sensor to determine rotational angles.
- Implemented algorithms to convert sensor readings into meaningful orientation information.




Figure 3.2: XDK Workbench

3. Integration of Dual Sensors:

- Integrated two XDK110 sensors S1 and S2 into the system: one for monitoring the collimator and another for environmental calibration.
- Designed and implemented a method to combine data from both sensors to accurately determine the orientation of the collimator.

4. Testing:

- A wooden replica of x ray collimator and its system was constructed for testing using wooden beams and angle connectors.



Figure 3.3 x-ray collimator replica for testing. [8]

5. Visualisation and Error Analysis:

- Created a Python-based visualization tool to display the orientation data in real-time, facilitating easier interpretation of results.
- Documented and analyzed the degree of error associated with the sensor measurements to assess the system's accuracy and reliability.




Figure 3.4: Collimator Rotation Principal [6]

4.2.3 Tools / Methods Used

1. Bosch XDK110 Sensor:

- Utilized for its capability to measure orientation. The sensor's built-in gyroscope and accelerometer were central to the project.

2. Programming Languages:

- **C/C++:** For initial code development and integration with the XDK110 sensor.
- **Python:** For developing the visualization tool and performing data analysis.

3. Development Environment:

- Bosch XDK workbench compatible with Bosch's XDK for sensor programming.
- Python libraries such as Matplotlib and NumPy for data visualization and analysis.

4. Data Analysis Techniques:

- Statistical methods to analyze the accuracy and error margins of sensor data.
- Algorithmic approaches for sensor fusion to combine data from multiple sensors.



Figure3.5: Tools Used.[7]

4.2.4 Project Goals

The primary goal of this project was to design and implement a robust system for detecting and visualizing the rotation angulation of a x-ray collimator. Collimators are critical components in medical imaging and radiation therapy, where precise alignment directly impacts the accuracy and effectiveness of diagnostic and therapeutic procedures.

To achieve this objective, the project utilized the Bosch XDK110 sensor, a versatile device known for its precision in motion detection. The specific Goals were to:

1. **Develop a Measurement System:** Create a system capable of accurately measuring the rotation angulation of the collimator using the Bosch XDK110 sensor. This includes calibrating the sensor to ensure precise data collection and integrating it into the collimator setup.
2. **Enhance Accuracy:** Implement additional sensors and data fusion techniques to improve the accuracy and reliability of the rotation measurements. This involves integrating multiple sensor inputs to cross-verify and enhance the precision of the system.
3. **Visualization and User Interface:** Design an intuitive visualization platform that effectively displays the rotation angulation data in real-time. The goal is to provide users with clear, actionable insights into the collimator's orientation, aiding in adjustments and alignment tasks.
4. **System Integration and Testing:** Ensure seamless integration of the measurement system with existing collimator setups and thoroughly test the system under various conditions to validate its performance and robustness.

4.2.5 Project Approach

1. Initial Setup and Familiarisation:

- Gain a comprehensive understanding of the Bosch XDK110 sensor.
- Set up the sensor in a controlled environment to assess its functionality.

2. Development of Orientation Visualisation:

- Write and test code to convert raw sensor data into angular orientation.
- Develop initial visualization tools to display orientation data.

3. Dual Sensor Integration:

- Implement integration of two sensors into the system.
- Develop and test algorithms for data fusion to improve accuracy.

4. Visualisation and Error Analysis:

- Create advanced visualization features in Python.
- Perform extensive testing to document and analyze sensor error margins.
- Generate comprehensive reports on system performance and accuracy.

4.2.6 Testing Stage



Figure 3.6: Test Setup of collimator angulation rotation.




Figure 3.7: Visualization on 3D plot.

4.2.6 Conclusion

The internship project on Collimator Rotation Angulation Detection was a successful endeavor, showcasing the effective application of the Bosch XDK110 sensor in a practical scenario. The project involved a comprehensive understanding of sensor technology, development of innovative software solutions, and integration of multiple sensors for enhanced accuracy.

4.2.7 Key Achievements include:

- **Effective Sensor Utilization:** Gained hands-on experience with the Bosch XDK110 sensor and utilized it for precise orientation measurement.
- **Functional Visualization Tools:** Developed a Python-based tool for real-time data visualization, significantly aiding in the interpretation of sensor readings.
- **Enhanced Accuracy:** Successfully integrated two sensors to improve the accuracy of orientation detection, with detailed documentation of error margins.

This project not only demonstrated technical proficiency in sensor integration and data visualization but also underscored the importance of accuracy in critical applications such as collimator alignment. The experience gained through this internship has provided valuable insights into both theoretical and practical aspects of sensor technology and data analysis.

5 Professional Insights

Internship Evaluation Involvement in various projects and activities throughout the internship period at ITT MP has provided an abundant of knowledge and insights in both technical and professional aspects.

5.1 Technical and Practical Knowledge

Technical requirements, in terms of knowledge in the topic and skillset for applications and software used, for each project vary. Through involvements in different types of projects such as Collision Avoidance System, product design and development, and prototype development for proof of concept and use-case testing, I am constantly challenged to learn and use different tools and technologies in a short span of time. With the continuous addition of technology, the ability to grasp a topic and implement it as part of your work in a timely manner is an important skill to have in a workplace. Prior to the internship, we were equipped with theoretical and practical knowledge that are fundamental as a mechatronics engineer. Although this is a good starting point, through this internship, I have learned and gained experience in new technologies and tools that are beneficial for my professional career as a mechatronics engineer. Finally, entrusted with several different project challenges me with both my time and project management skills. The ability to prioritize and portion the amount of effort and work given to each project or task to ensure that every project is moving forward and ready to be presented by the next update meeting is scheduled. In addition, communication within the project team is key to make certain that each task is handled and delegated accordingly.

5.2 Non-technical outcomes

At Siemens Healthineers, preparing documentation adheres to established standards, while collaboration within an internationally diverse team is crucial. Understanding workplace conduct and safety regulations is essential for maintaining a productive office environment. Efficiently managing multiple tasks and meeting deadlines are key responsibilities, alongside

participating in tests and learning opportunities through ITT Biomedical Technologies discussions. Gaining insight into corporate working methods and the ability to learn and collaborate on various projects with teammates are also vital components of the role.

5.3 Conclusion

When I think about my time at Siemens Healthineers, I remember the friendly atmosphere and how everyone worked together. My colleagues were welcoming, and we worked as a team to make our projects successful. One important thing I learned was how important communication is. I talked a lot about our projects, divided tasks, and had meetings with others. This helped us make progress and taught me that good communication means both speaking and listening. This internship also showed me that learning never stops. Whether I watched my colleagues in meetings or learned to use new tools, each experience made me better. I realized that being open to new things and adapting to change is important for growing. Reflecting on my internship journey at Siemens Healthineers, I am filled with a deep sense of accomplishment and gratitude. This chapter of my professional life has been a remarkable blend of learning, growth, and meaningful experiences. As I conclude this phase, I want to express my sincere appreciation for the invaluable lessons and opportunities I've been fortunate to receive. Moreover, my abilities in programming, prototyping, managing projects and tasks, and collaborating within teams have significantly expanded. Additionally, I had the chance to put into practice the knowledge I gained in school, particularly in areas such as hardware Engineering, programming, design and printing. Working closely with colleagues, engaging in open discussions, and adapting to changing scenarios have not only enhanced my teamwork skills but have also instilled in me the importance of a growth mindset.

6 Certificates at ITT MP

7 Internship Certificate

8 Bibliography

- [1] “Siemens Healthineers | Corporate Home.” Accessed: Aug. 06, 2024. [Online]. Available: <https://www.siemens-healthineers.com/>
- [2] “Siemens Healthineers invests €60 million in Kemnath production location.” Accessed: Aug. 06, 2024. [Online]. Available: <https://www.siemens-healthineers.com/press/releases/kemnathinvest>
- [3] “Nexaris Angio-CT.” Accessed: Aug. 06, 2024. [Online]. Available: <https://www.siemens-healthineers.com/angio/nexaris-therapy-suites/nexaris-angio-ct>
- [4] “operating_instructions_nanoscan3_i_o_en_im0087137.pdf.” Accessed: Aug. 06, 2024. [Online]. Available: https://cdn.sick.com/media/docs/7/37/137/operating_instructions_nanoscan3_i_o_en_im0087137.pdf
- [5] “Safety Designer - Safety Designer | SICK.” Accessed: Aug. 06, 2024. [Online]. Available: <https://www.sick.com/de/de/catalog/produkte/safety/sicherheitssteuerungen/safety-designer/safety-designer/p/p674218>
- [6] “MULTIX Impact C mit myExam Companion.” Accessed: Aug. 06, 2024. [Online]. Available: <https://www.siemens-healthineers.com/de/radiography/digital-x-ray/multix-impact-c>
- [7] “XDK_Guide_General_Information.pdf.” Accessed: Aug. 06, 2024. [Online]. Available: https://developerportal.xdk.io/downloads/XDK_Guide_General_Information.pdf
- [8] “X-ray Collimator.” Accessed: Aug. 06, 2024. [Online]. Available: <https://www.oem-products.siemens-healthineers.com/x-ray-collimator>